


## ● PathoProof<sup>®</sup> Mastitis PCR Assay



Martina Kahila, Product Manager PathoProof products

# The World Leader in Serving Science

## ThermoFisher S C I E N T I F I C

### Who We Are

**Leading provider of analytical instruments, equipment, reagents and consumables, software and services for research, analysis, discovery and diagnostics**

### Strengths

<b>Scale</b>	<ul style="list-style-type: none"><li>• \$11 billion revenues</li><li>• 34,000 employees</li><li>• 350,000 customers</li><li>• 150 countries served</li></ul>
<b>Capabilities</b>	<ul style="list-style-type: none"><li>• Complete portfolio</li><li>• World-class technologies</li><li>• Commercial reach</li></ul>
<b>Experience</b>	<ul style="list-style-type: none"><li>• 150 years of combined experience</li></ul>
<b>Brand Equity</b>	<ul style="list-style-type: none"><li>• Pre-eminent brands</li></ul>

***Enabling customers to make the world healthier, cleaner and safer***

\* Figures from Thermo Fisher Scientific Corporate Profile



## In research and analysis

- Advanced analytical instrument technologies
- Solutions to improve laboratory workflows



## In healthcare and clinical laboratories

- Diagnostic kits, reagents and consumables
- Clinical trial support and total supply



## In manufacturing and the field

- Safety-related products and personal protection
- Environmental and process instruments

***World-class analytical technologies and the most complete portfolio of laboratory products and services***

● Microbiology Division



# About the Microbiology Division:

- Specialist business focused on Microbiology for diagnosis of infectious disease and detection of bacterial contamination in industrial applications.
- Legacy brands of Oxoid and Remel, world renowned in culture media and diagnostic assays, combining with Thermo Scientific brand which represents technology and innovation.
- Long-standing collaboration with DuPont Qualicon to commercialise the BAX range of PCR products for food safety in Europe, Australia and NZ.
- Direct sales and service in 24 countries
- Focused R&D with an expanded Finnzymes team to continue to grow the product range




*Combining expertise in the science of microbiology with leading edge technology to continue to develop better ways to detect and identify micro-organisms`*

# Mastitis

---

- Responsible for huge losses each year to the dairy industry due to many reasons, for example reduced milk production and discarded milk unfit for human consumption
- Economical losses ~200 USD/cow/year. Total losses in the USA alone over 1,7 billion dollars per year
- Bacterial culturing is still regarded as the golden standard in mastitis diagnostics
- However, there are many problems related to culturing, for example:
  - Slow
  - Unspecific
  - Sensitive to errors
  - Can not be used with preserved milk samples

# PCR=Polymerase Chain Reaction



10 cycles:  $2^{10} \sim 1000$  fold amplification

20 cycles:  $2^{20} \sim 1000\ 000$  fold amplification


30 cycles:  $2^{30} \sim 1000\ 000\ 000$  fold amplification

# Real-time PCR

- PathoProof is based on real-time PCR, which is the new golden standard in for example human clinical diagnostics and food diagnostics
- When using real-time PCR, the amplification and the detection of DNA is done simultaneously
- Real-time PCR is one-step, single tube detection, and the sample tube stays closed during amplification, minimizing contamination risks
- Real-time PCR is a quantitative method, which means that the initial amount of each bacteria can be detected

Fluorescence curves represent the bacterial amounts in the sample

C(t) value = the PCR-cycle where the fluorescence grows over a certain threshold level



# PathoProof Mastitis PCR Assay

- Fast
- Sensitive
- Accurate
- Easy to use
- Can be used with:
  - Bulk tank milk samples
  - Individual quarter milk samples
  - DHI samples



# Speed, accuracy and sensitivity

- 3-4 hours from sample to result
- 100 % specificity – the test detects the bacterial DNA and interpretation of results is done by a software designed for PathoProof
- No 'mixed growth' results. All bacteria are detected simultaneously. The most prevalent pathogen is reported if it accounts for >90% of the total bacterial quantity
- 40-50 % of samples yielding no growth in conventional culturing provide positive results with the PathoProof Mastitis PCR Assay
- Equally sensitive for viable and growth inhibited bacteria



# PathoProof vs. culturing



Benefits of PathoProof over bacterial culturing in mastitis diagnostics are tremendous	
Bacterial culturing	PathoProof Mastitis PCR Assay
48-72 hours	4-5 hours for all pathogens, incl. <i>Mycoplasma</i>
25-40% of no-growth cases	Detects growth-inhibited and dead bacteria -significantly less (>40% less) 'no growth' results
User experience affects reliability	Objective results for all bacteria -after product training high accuracy can be reached in all laboratories
Mainly qualitative	Qualitative and quantitative -calculates bacterial proportions in 'mixed growth' cases
Unreliable during antibiotic treatment	Not affected by antimicrobials -can be used for following and studying treatment efficacy -reliable for repeated testing even if antibiotics have been/are being used
Prone to error due to bacterial growth during sample transportation	Can be used with bronopol-preserved milk (elimination of error)
Cannot be used for testing of DHI samples	Can be used for testing of DHI samples

# Easy to use

- Laboratories with no previous microbiology or molecular biology experience use the assay successfully and get equally accurate results
- Kits include all DNA extraction and PCR reagents
- The software automatically analyze the results and create results reports
- We provide a 3-4 days training to use the kits and the instruments. The training is free of charge in combination to the real-time PCR instrument purchase








**1. DNA-extraction from fresh or preserved milk samples**




**2. Real-time PCR amplification of bacterial DNA**

**3. Simultaneous detection and quantification of the 11 most important mastitis species/groups and the most relevant antibiotic resistance gene (bla<sub>Z</sub>)**



**4. Analysis of results**



# PathoProof kit configurations


## PathoProof Mastitis Complete-12

*Staphylococcus aureus*  
*Coagulase negative staphylococci*  
*Streptococcus agalactiae*  
*Streptococcus dysgalactiae*  
*Streptococcus uberis*  
*Escherichia coli*  
*Corynebacterium bovis*  
*Enterococcus faecalis and faecium*  
*Klebsiella pneumoniae and oxytoca*  
*Serratia marcescens*  
*Arcanobacterium pyogenes and Peptostreptococcus indolicus*  
Staphylococcal beta-lactamase gene

⇒ **These bacteria are responsible for >95 % of all clinical and subclinical mastitis cases**


## PathoProof Mastitis Major-3

*Mycoplasma bovis*  
*Staphylococcus aureus*  
*Streptococcus agalactiae*



**PathoProof™  
Mastitis PCR Assay**

Revolutionary  
fast and reliable  
identification of  
mastitis causing  
bacteria from  
bovine milk



**FINNZYMES  
DIAGNOSTICS**

# PathoProof Software

- Automatically analyzes the PCR results
- Creates reports, which show the amount of each bacteria separately
- If one bacteria is present in over 90 % of the total bacterial load, then that is also reported

Category	Explanation
-	Bacterial DNA not detected
+	Bacterial DNA detected in small amounts
++	Bacterial DNA detected in intermediate amounts
+++	Bacterial DNA detected in large amounts

# Field study: PCR vs. bacterial culturing

Koskinen et. al., *J.Dairy Sci.*, In press

- First comprehensive field comparison about PCR and culturing from mastitic milk samples
- Included 1000 milk samples from cows with clinical or subclinical mastitis or from clinically healthy quarters with low SCC
- 77 % of the clinical samples were positive in culturing and 89 % in PCR
- "Healthy samples" provided a negative result in most cases (4,3 % of samples positive in culturing and 7,5 % in PCR), indicating that PCR is not "too sensitive"
- Benefits of PCR were clear: speed, automated interpretation of results as well as increased sensitivity



© 2014 ThermoFisher Scientific

0000000000000000

# PathoProof customers

---

- The PathoProof system has been installed into about 30 laboratories in 16 countries
- Currently about 60 % of the laboratories are DHI laboratories and 40 % animal health laboratories
- During the first 9 months in 2010, about 90 % of all installations have been to DHI laboratories
- The highest throughput is currently in Finland at Valio laboratories, where about 100 000 samples are analyzed per year
- Our service laboratory in Finland analyzes about 15 000 samples per year
- We organized the first PathoProof User's Day on September 1st, 2010 in Finland. This will be a yearly event in the future.

# PathoProof feedback from laboratories and their customers

- Convenience is very important, many laboratories use the DHI samples to analyze mastitis
- Different sample types are being used: individual quarter milk samples, string samples, combined quarter samples, DHI samples, bulk tank milk samples..
- More expensive, but:
  - Lower cost of sample transportation
  - Less resampling (no growth results)
  - Faster results (the same day)
  - Integration with the DHI results
  - SCC results can be utilized efficiently
- In Finland, Valio made a survey about the milk producers interest in PathoProof testing, 98 % answered that they want their samples tested with PCR



# PathoProof kit formats

- Large kits (F-870L and F-914L)
  - Reagents for 4 x 96 samples
  - DNA extraction is done on 96-well plates
- Small kits (F-870S and F-914S)
  - Reagents for 50 samples
  - DNA extraction is done in single tubes



---

Thank you for your attention!

Questions?